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Report, Bookmark or events, Averages, Write Screen

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Worst case, Periodic average, peak exercise

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Protocol name, Heart rate (HR)/% of target HR, ST level and ST slope, BP,

Stage and Stage time, METs, Double Product, One average complex, Running

trends of Me Treadmill speed and grade or watts ergometer

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USAF/SAM 3.3, Low Ramp, Medium Ramp, High Ramp, Time Ramp, METs Ramp,

Astrand(ergometer), Pharmacological

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
# Q-Stress

## PRIMARY BENEFITS:

- **Logical user interface:** The system's left-to-right navigation and intuitive controls guide users through procedures with ease.
- **Improved trace quality:** Our innovative Source Consistency Filter (SCF) minimizes noise and baseline artifact while preserving a diagnostic-quality ECG signal.
- **Wireless freedom:** The Wireless Acquisition Module (WAM™) digitally transmits ECG data, providing comfort and freedom of movement for the patient.
- **12-Lead ECG interpretation:** Best-in-class VERITAS™ algorithm provides 12-Lead ECG interpretation for adults, adolescents and children.

3.1.14

## A UNIQUE APPROACH TO CLINICAL DATA

Q-Stress®, brought to you by Mortara, is designed to provide comprehensive clinical information to meet the ever-changing demands of today's diverse patient population, from routine exams to complex cardiac evaluations.

- Quickly access continuous ST trends and profile, ST level and slope for all 12 average complexes, and auto comparison with reference complexes.
- Enhanced full disclosure mode provides a beat-by-beat review during the exam, as well as post-test, allowing easy retrospective addition of events and ECG printing.

## A SOLUTION DRIVEN DESIGN TO ENHANCE YOUR WORKFLOW

Designed by clinicians for clinicians, the powerful graphical user interface paired with our 24-inch widescreen color touch display works with you - the way you work.

- Begin by importing your patients' demographics. This fast and simple step enables you to eliminate manual data entry, streamline your work, and minimize documentation and billing errors.
- Append significant data points as needed to the final report.
- Configured with our TM55 or TM65 treadmill, Q-Stress will ensure you're fully equipped for top performance. Our medical treadmills have been the industry's gold standard for more than three decades and are among the best selling in the world.

## A CONNECTIVITY OPTION THAT'S RIGHT FOR YOU

From simple connectivity to complete interoperability, Q-Stress offers a variety of integration solutions.

- Communication choices include PDF, XML, HL7, and DICOM®.
- Combine stress data with other systems that use non-proprietary formats, including your echo and nuclear imaging systems.
- Share final reports, access shared printers, and archive information to your EHR.

12-Lead ECG interpretation: Best-in-class VERITAS™ algorithm provides 12-Lead ECG interpretation for adults, adolescents and children.

QUINTON

# QUINTON

## MORTARA INSTRUMENT, INC.

10000 W. 16th Ave.  
Minneapolis, MN 55426  
Tel: 612-835-1000  
Fax: 612-835-1001  
E-mail: sales@mortara.com

## MORTARA INSTRUMENT EUROPE S.R.L.

Via S. Giovanni Battista 1  
20121 Milano, Italy  
Tel: +39 02 574011  
Fax: +39 02 574012  
E-mail: sales@mortara.com

## MORTARA INSTRUMENT AUSTRALIA

PO Box 100  
Bentley, Victoria 3084  
Australia  
Tel: +61 3 9401 1000  
Fax: +61 3 9401 1001  
E-mail: sales@mortara.com

## MORTARA INSTRUMENT UK LTD

Unit 10, The Gateway  
St. Mary's, Basingstoke, Hampshire  
RG24 0AP, UK  
Tel: +44 1256 334444  
Fax: +44 1256 334445  
E-mail: sales@mortara.com

www.mortara.com

ISO 13485 CERTIFIED

## WARRANTY & SERVICE

Mortara Instrument is committed to the highest level of customer support. Please contact us for the program which best suits your needs.

\*Specifications subject to change without notice

Burdick®, Mortara®, Q-Stress®, Quinton®, VERITAS™ and WAM™ are trademarks or registered trademarks of Mortara Instrument, Inc.

Microsoft® and Windows® are registered trademarks of Microsoft Corp.

SunTech™ and Tango™ are registered trademarks of SunTech Medical, Inc.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information

2016-01

# Q-STRESS®

FEATURE	SPECIFICATION*
<b>OPERATING SYSTEM</b> 3.1.1	Microsoft® Windows® 7 Professional 64-bit
<b>POWER REQUIREMENTS</b>	100/120 VAC 50/60 Hz 2.5 A nominal; 200/240 VAC 50/60 Hz 1.3 A nominal
<b>PROTOCOLS</b> 3.1.13	<b>STANDARD</b> – Bruce, Modified Bruce, Modified Balke, Naughton, USAF/SAM 2.0, USAF/SAM 3.3, Ramped-low, Ramped-medium, Ramped-high, Time-Ramp, and HES Ramp, Astrand (ergometer), pharmacological and more; <b>CUSTOM</b> – Up to 100 custom protocols supported
<b>FRONT END MODULE</b>	AM12Q with AHA 10-electrode-pinch or snap connectors AM12Q with IEC 10-electrode-pinch or snap connectors WAM™ for wireless ECG acquisition with AHA 10-electrode-pinch connectors WAM for wireless ECG acquisition with IEC 10-electrode-pinch connectors
<b>ECG</b> 3.1.4 3.1.12.2	<b>LEAD GROUPS</b> – Standard (Mason-Likar), Cabrera 12-lead <b>DISPLAY FORMAT</b> – 3, 6, 6x2, or 12x1 channels <b>PERFORMANCE STANDARD</b> – AAMI-EC11: gain accuracy, frequency response, CMRR, system noise, dynamic range, and input impedance
<b>DISPLAY AND ANALYSIS FILTERS</b>	Source Consistency Filter (SCF), muscle artifact, baseline wander, 40 Hz low pass, line frequency
<b>REPORTS</b> 3.1.11, 3.1.11.1 3.1.11.2	<b>IN-TEST</b> – 12-lead, average beat report, 1-page write screen, continuous rhythm, arrhythmia events <b>FINAL</b> – Summary, tabular, worst case, average beats by stage or by minute, in-test ECG reports, trend graphs and peak exercise. Reports can be viewed, printed, saved in PDF format, emailed or exported in PDF, XML and DICOM. <b>CUSTOM</b> – Unlimited customizable reports
<b>OUTPUTS</b> 3.1.8	12-lead pulse and analog signals for QRS detection; synchronization using any ECG lead
<b>EXTERNAL DEVICE INTERFACES</b>	<b>TREADMILLS</b> – Quinton TM 55/65, Trackmaster TMX428 and TMX425 <b>ERGOMETERS</b> – Compatible with Ergoline and Lode communication protocols <b>BLOOD PRESSURE MONITORS</b> – SunTech® Tango+ and Tango M2
<b>PRINTERS</b>	<b>THERMAL</b> – High-resolution Z200+ with automatic feed and continuous printing capability (requires additional network port) <b>LASER</b> – Optional compact HP M401n laser printer with high quality print (recommended) <b>PAPER</b> – 8.5 x 11" (letter) or 210 x 300 mm (A4)
<b>NETWORK INTERFACE</b>	Microsoft-compatible networking for storage, distribution and email; 100 Mbps connection or better is required
<b>EXPORT/COMMUNICATION PROTOCOL/FORMAT</b> 3.1.19	TCP/IP, PDF, XML, HL7, DICOM
<b>REMOTE TECHNICAL SERVICE INTERFACE</b>	Remote service increases system uptime availability and decreases the user's requirement to assist the troubleshooting process. Available via network or desktop connections (only in the US).
<b>DISPLAY</b>	24" (609 mm) flat panel LCD; 24" touchscreen display is optional
<b>CART DIMENSIONS AND WEIGHT</b>	36 x 30 x 20" (914 x 762 x 508 mm); 150 lb (68 kg); (weight varies based on cart options)
<b>WARRANTY</b>	13 months parts and labor
<b>OPTIONS</b>	Networking, Resting ECG Analysis, ST-segment Re-Analysis, TTL Trigger, WAM

## ORDERING INFORMATION – ACCESSORIES

9922-020-50  
9922-017-52  
9100-026-60  
9100-026-61

SunTech Tango M2 BP Monitor Kit for Q-Stress  
SunTech Tango and Tango M2 SpO2 option  
Thermal Paper, Z200+, Z-fold, full grid, 8.5 in x 11 in, 12 packs/case  
Thermal Paper, Z200+, Z-fold, full grid, A4 format, 210mm x 300mm, 12 packs/case

BURDICK Mortara QUINTON

Ref: 9513-192-01-ENG Rev 4

### Waveform Control and Using Filters

Left-click anywhere in the real-time ECG to open a window allowing you to set displayed ECG leads, filters, display gain, and display speed.

Filters listed below may be enabled or disabled at any time during the exam:

3.1.2

- Source consistency filter (SCF)
- 40 Hz filter
- AC filter

**Waveform Control**

**Lead layout**

- ☐ 3-Lead
- ☐ 6-Lead
- ☒ 6x2-Lead
- ☐ 12-Lead

**Filters**

- ☒ AC Filter
- ☐ 40 Hz filter
- ☒ SCF

**Gain**

- ☐ 2.5 mm/mV
- ☐ 5 mm/mV
- ☒ 10 mm/mV
- ☐ 20 mm/mV
- ☐ 40 mm/mV

**Display speed**

- ☒ 25 mm/s
- ☐ 50 mm/s

3.1.3

### Source Consistency Filter (SCF)

Mortara Instrument's patented Source Consistency Filter (SCF) is an exclusive feature used for reducing noise associated with stress testing. Using the morphology learned during the Pre-exercise phase or during a relearn operation, the SCF distinguishes between noise and the actual signal in each of the 12 leads. This filtering reduces muscle tremor noise, low-frequency and high-frequency noise, and baseline artifact while preserving diagnostic-quality waveforms.

The default state of the filter (enabled or disabled) is determined by the modality settings. When the filter is on, SCF appears in the lower border of the real-time ECG display. This setting can be changed at any time during a stress exam.

### AC Filter

The AC filter removes line frequency noise by removing frequencies in a narrow band around 60 Hz (domestic) or 50 Hz (international). The default state of the filter (enabled or disabled) is determined by the selected profile. When the filter is on, 60 Hz or 50 Hz appears in the lower border of the real-time ECG display. This setting can be changed at any time during a stress exam.

### 40 Hz Filter

The 40 Hz filter is a plot filter, meaning that it affects only the plotted/printed information similar to a 40 Hz filter on an electrocardiograph. The default state of the filter (enabled or disabled) is determined by the selected profile. When the filter is on, 40 Hz appears in the lower border of the real-time ECG display. This setting can be changed at any time during a stress test.

**WARNING:** When the 40 Hz filter is used, the frequency response requirement for diagnostic ECG equipment cannot be met. The 40 Hz filter significantly reduces high-frequency components of the ECG and pacemaker spike amplitudes, and is recommended only if high-frequency noise cannot be reduced by proper procedures.

## 24. Q-STRESS PHYSICIAN'S GUIDE

---

### Q-Stress Signal Analysis

This appendix of the user manual is intended for qualified health care professionals who want to learn more about the Q-Stress cardiac stress testing system methods for analysis. It describes various aspects of the Mortara VERITAS™ algorithm for cardiac stress testing analysis.

The Q-Stress analysis process consists of the following steps and components:

1. Acquiring the raw digital ECG data and then processing and storing the digital ECG data.
2. Beat detection and creation of the dominant QRS template
3. Determination of Normal, ventricular, and paced beat types
4. Resting ECG analysis during the Pre-Exercise phase
5. Identification of dominant QRS change in rhythm
6. ST segment analysis of all 12-leads
- 3.1.15 7. **Arrhythmia detection**
8. Cardiac stress testing calculations and algorithms
9. Creation of summary data to provide results reporting

### Data Acquisition

3.1.5

The Q-Stress patient cable acquires the ECG signal that is digitized at a sampling rate of 40,000 samples per second per channel by the Q-Stress front end (AM12Q). The AM12Q also provides a mode of operation for measuring the patient electrode impedance. The ECG data is filtered and stored by the Q-Stress system at a sampling rate of 500 samples per second per lead, with a resolution of 2.5 microvolts per LSB. The frequency bandwidth of the ECG signal meets AAMI standard EC11:1991 for diagnostic quality ECG data.

### Full Disclosure Display and Storage

Full disclosure of historic ECG data is navigable through use of the Context View from the start of the Pre-exercise phase through Recovery end. A miniature single lead is user defined and can be changed to any of the 12-leads at any time. The Context View area highlight stored ECG events and allows addition of new ECG events as well as deletion of stored events that are not desired.

Post exam review supports full disclosure data navigation beat-by-beat as well as quick display of any event from Pre-exercise to Recovery end by a click on ECG events that occurred during the exam, or by a selection of any point in the trends.

### Final Report


Cardiac stress testing results can be reported and exported in PDF, XML, and DICOM formats. The final report pages are numbered sequentially with examples explained in this guide.

**Lead Impedance Check 3.1.6**

Before moving to the pre-exercise phase, electrodes may be checked to verify a satisfactory connection. The Q-Stress Impedance Check is recommended during the Observation phase.

Select the **Impedance Check** button to have the program check the quality of each electrode connection. A window will open showing your electrode connections.

*NOTE: Impedance Check is only available when the Q-Stress patient cable is selected as the front end.*

Color	Connection Quality	Impedance (kΩ)
Green	Good	Less than 75
Yellow	Fair	75 – 250
Red	Poor	Greater than 250
White	Open (disconnected)	--

Re-prepare the site(s) and replace the electrode(s) until you have confirmed good connection, then select the Exit button.

*NOTE: There is a 5-second delay where no waveform is present on the observation screen after exiting the Impedance Check display and returning to real-time ECG display.*

**Demonstration Mode** 3.1.7


The Q-Stress includes a demonstration mode that provides a way to present features and train clinicians on system operation without the need for live patient connection.

The demonstration mode is enabled when a patient name of Demo and a number or multiple numbers is entered in the Last Name field (e.g. Demo1 or Demo2 or Demo123 and so on). The word "Demonstration" is watermarked within the ECG waveform to differentiate the display from live patient ECG waveform.

**NOTE:** The *D* must be upper-case and *emo* must be lower-case or the demonstration mode will not be enabled.

The user interface and display in the demonstration mode is no different than the user interface and display in the live mode, with a few exceptions:

- BP readings with a configured Tango BP monitor are not initiated. Demo BP values are periodically displayed and updated during an exam.
- Configured exercise equipment (treadmill or ergometer) is not controlled when in the demonstration mode.

# Features & Technical Specifications


Feature	OptiPlex XE3 Technical Specifications <sup>1</sup>
Processors <sup>1</sup>	Intel Celeron™ G4900 (2 Cores/2MB/2T/3.1GHz/65W); supports Windows 10/Linux Intel Pentium™ Gold G5400 (2 Cores/4MB/4T/3.7GHz/65W); supports Windows 10/Linux 3.1.9 Intel® Core™ i3-9100 (4 Cores/8MB/4T/3.6GHz/65W); supports Windows 10/Linux Intel® Core™ i5-8500 (6 Cores/9MB/6T/up to 4.1GHz/65W); supports Windows 10/Linux Intel® Core™ i7-8700 (6 Cores/12MB/12T/up to 4.6GHz/65W); supports Windows 10/Linux
Chipset	Intel® Q370 Chipset
Operating System Options <sup>1</sup>	Microsoft® Windows 10 Home (64-bit) 3.1.1 Microsoft® Windows 10 Pro (64-bit) Ubuntu® 16.04 LTS (64-bit)
Graphics Options <sup>2</sup>	Integrated Intel® HD Graphics 630 AMD Radeon™ R5 430, 2GB AMD Radeon™ RX 550, 4GB NVIDIA GeForce® GT 730, 2GB NVIDIA GeForce® GTX 1050, 4GB (Tower only) Dual AMD Radeon™ R5 430, 2GB Dual AMD Radeon™ RX 550, 4GB (Tower only)
Memory	4 DIMM slots; Non-ECC dual-channel 2666Hz DDR4 SDRAM (Memory performance on Intel® Core i3 processor and below will be at 2400MHz). Max memory is 64GB 3.1.9 4GB (1x4GB) 2666MHz DDR4 Memory 8GB (2x4GB) 2666MHz DDR4 Memory 8GB (1x8GB) 2666MHz DDR4 Memory 16GB (2x8GB) 2666MHz DDR4 Memory 16GB (1x16GB) 2666MHz DDR4 Memory 32GB (4x8GB) 2666MHz DDR4 Memory 32GB (2x16GB) 2666MHz DDR4 Memory 64GB (4x16GB) 2666MHz DDR4 Memory
Networking <sup>1,4</sup>	Integrated Intel I219-LM Ethernet LAN 10/100/1000 Qualcomm® QCA61x4A Dual-band 2x2 802.11ac Wireless with MU-MIMO + Bluetooth 4.2 Intel® Wireless-AC 9560, Dual-band 2x2 802.11ac Wi-Fi with MU-MIMO + Bluetooth 5
I/O Ports	Tower/SFF 10 External USB: 1 x USB Type-C 3.1 Gen 2; 5 x 3.1 Gen 1 Type-A (1 front, 4 rear) and 4 x 2.0 (2 Type-A front (one with PowerShare), 2 Type A rear (both SmartPower On)) 1 RJ-45 2-Displayport 1-Serial 2-PS/2 1 Optional 3rd Video Port (VGA/DP/HDMI 2.0b/USB Type-C alt-mode) Optional Powered USB add-in card, Optional Powered Serial add-in card 1 Universal Audio Jack 1 Line-Out
Removable Media Options	Optional 8x DVD-ROM 9.5mm ODD Optional 8x DVD+/-RW 9.5mm ODD Optional SD Media Card Reader
Storage Options <sup>1</sup>	2.5 inch 500GB 5400rpm SATA Hard Disk Drive 2.5 inch 500GB 7200rpm SATA Hard Disk Drive 2.5 inch 500GB 5400rpm SATA Solid State Hybrid Drive w/ 8GB Flash 2.5 inch 500GB 7200rpm FIPS Self Encrypting Opal 2.0 Hard Disk Drive 2.5 inch 1TB 7200rpm SATA Hard Disk Drive 2.5 inch 1TB 5400rpm SATA Solid State Hybrid Drive w/ 8GB Flash 2.5 inch 2TB 5400rpm SATA Hard Disk Drive 3.5 inch 500GB 7200rpm SATA Hard Disk Drive 3.1.9 3.5 inch 1TB 7200rpm SATA Hard Disk Drive 3.5 inch 2TB 7200rpm SATA Hard Disk Drive

## Q-Stress Specifications


Feature	Workstation Minimum Specification*
Computer Model	HP rp5800 CPU
Processor	Intel Core i3
Graphics	1920 x 1080 or 1920 x 1200
RAM	4 GB
Operating System	Microsoft Windows 7 Professional 32-bit or 64-bit
Hard Drive Capacity	500 GB
Archive	Network or external USB drive
Input Devices	Standard keyboard and scroll mouse
Software Installation	CD-ROM
Network	100 Mbps connection or better
Front-end ECG Devices	AM12Q with Patient Cable connector Wireless Acquisition Module (WAM) Trigger Module for analog and TTL signal output to external devices
Printing Devices	HP M401n Laser printer (recommended) Z200+ thermal printer (requires additional network port)
USB Ports	2 free USB 2.0 ports
Serial Ports	2 serial ports for serial interfaced equipment
Audio	Required for NIBP and Pharmacological notification
<b>Isolation Transformer - Required When Workstation is Used for Stress Testing</b>	
Isolation Transformer Requirement 3.1.10	Known Agency Mark (KAM) Meets IEC 60601-1 requirements Protective Earth conductor for all connected equipment Z200+ only configuration: 300 watts LaserJet printer configuration: 1,000 watts
Feature	Server Minimum Specification*
Processor	Performance equivalent to an Intel Xeon class, Quad-core with hyperthreading
Graphics	1024 x 768
RAM	4 GB
Operating System	Microsoft Windows 2008 or 2012 server R2, 64-bit
System Disk	100 GB for OS and product installation (RAID recommended for data redundancy)
Data Disks	550 GB hard drive space available HD controller with 128 MB read/write cache (RAID recommended for data redundancy)
Archive	Network or external USB drive
Software Installation	CD-ROM
Network	100 Mbps connection or better
Input Devices	Standard keyboard and mouse

\* Specifications subject to change without notice.

Designed with essential features  
for home and business

Maximum tilt 5.5° forward or 22° backward.

#### Upgrade to P2314H<sup>5</sup>

In-plane switching panel technology.  
DisplayPort connectivity

#### What's in the box?

Monitor, stand, VESA screw cover, power cable, drivers, documentation media, Quick Set Up Guide and Safety Information, VGA cable, DVI-D cable.

Dell recommends that customers dispose used computer hardware, including monitors, in an environmentally sound manner. Potential methods include reuse of parts or whole products and recycling of product, components and/or materials. For more information, please visit [http://dell.com/recycling\\_programs](http://dell.com/recycling_programs) and [www.dell.com/environment](http://www.dell.com/environment)

## Dell E2314H (58.42 cm; 23" VIS) 3.1.12

### Display

Model number	E2314H
Viewable image size (diagonal) <sup>1</sup>	58.42 cm (23 inches) 3.1.12.1
Active display area	
Horizontal	509.8 mm (20.07 inches)
Vertical	286.7 mm (11.29 inches)
Maximum resolution	Full HD 1920 x 1080 at 60 Hz
Aspect ratio	16:9
Pixel pitch	0.255 mm
Brightness (typical)	250 cd/m <sup>2</sup>
Color gamut (typical)	83% <sup>2</sup>
Color depth	16.7 million colors
Contrast ratio (typical)	1000:1
Viewing angle (typical) (vertical/horizontal)	160°/ 170°
Response time (typical)	5 ms (black to white)
Panel technology	TN (active matrix - TFT LCD) 3.1.12.1
Backlight	LED

### Connectivity

Connectors	VGA, DVI-D (HDCP)
Dell soundbar <sup>4</sup> (optional)	ACS11
Remote asset management	Yes, via Dell Display Manager

### Design features

Stand	Tilt only (Typical: 5° forward or 21° backward; Maximum: 5.5° forward or 22° backward)
VESA mounting support (wall mount kit sold separately)	100 mm x 100 mm
Security	Security lock slot and stand lock (security lock not included)

### Power

AC input voltage/frequency/current	100 to 240 VAC/50 or 60 Hz ± 3 Hz/1.5 A (typical)
Power consumption	22W (typical)
Power consumption standby/sleep mode	<0.5 W

### Dimensions (with stand)

Height	405.4 mm (15.96 inches)
Width	546.2 mm (21.50 inches)
Depth	165.5 mm (6.52 inches)

### Weight

Weight (panel only - no stand)	3.13 kg (6.89 lb)
Weight (with stand)	4.22 kg (9.28 lb)
Shipping weight	5.50 kg (12.10 lb)

### Standard service plan

• 3 years Advanced Exchange Service<sup>3</sup> & Limited Hardware Warranty<sup>4</sup>

### Environmental compliance

ENERGY STAR<sup>®</sup>, TCO Certified Displays, EPEAT<sup>®</sup> Gold<sup>1\*</sup>, China Energy Label, WEEE, ErP (EuP) Standards, Korea E-Standby

<sup>1</sup> Registered in US and Canada only

<sup>2</sup> EPEAT Gold registered in the US. EPEAT registration varies by country. See [www.epeat.net](http://www.epeat.net) for registration status by country

<sup>3</sup> For a copy of the Limited Hardware Warranty, write Dell USA LP, Attn: Warranties, One Dell Way, Round Rock, TX 78682 or see [dell.com/warranty](http://dell.com/warranty).

<sup>4</sup> Advanced Exchange Service: Replacement part/unit dispatched, if needed, following completion of phone/online diagnosis. Fee charged for failure to return defective unit. Availability varies. Other conditions apply.

<sup>5</sup> Color gamut (typical) is based on CIE1976 (83%) and CIE1931 (72%) test standards.

<sup>6</sup> Only compatible with ACS11 soundbar. Launch date for soundbar varies across countries. Please check with your country sales representatives for more information.

<sup>7</sup> Please check with your sales representative for availability of upgrade product.

Based on Ad #A1300196B [dell.com/monitors](http://dell.com/monitors)

Product availability varies by country. Please contact your Dell representative for more information.

- **12-LEAD IN 6x2-FORMAT WITH ZOOMED ST ANALYSIS WITHOUT TRENDS**  
Twelve channels consisting of 6 seconds of twelve leads
- **12-LEAD IN 6x2-FORMAT WITHOUT ZOOMED ST ANALYSIS AND WITHOUT TRENDS**  
Twelve channels consisting of 8 seconds of twelve leads
- **12-LEAD IN 6x2-FORMAT WITHOUT ZOOMED ST ANALYSIS AND WITH TRENDS**  
Twelve channels consisting of 6 seconds of twelve leads
- **12-LEAD IN 12x1 FORMAT WITH ZOOMED ST ANALYSIS AND TRENDS**  
Twelve channels consisting of 8 seconds of twelve leads
- **12-LEAD IN 12x1 FORMAT WITH ZOOMED ST ANALYSIS WITHOUT TRENDS**  
Twelve channels consisting of 12 seconds of twelve leads
- **12-LEAD IN 12x1 FORMAT WITHOUT ZOOMED ST ANALYSIS AND WITHOUT TRENDS**  
Twelve channels consisting of 15 seconds of twelve leads
- **12-LEAD IN 12x1 FORMAT WITHOUT ZOOMED ST ANALYSIS AND WITH TRENDS**  
Twelve channels consisting of 12 seconds of twelve leads

The Context View can be enabled to show a single lead from start of the pre-exercise phase through the recovery phase indicating automatic and manually stored ECGs. This view also allows addition of historic ECG events and deletion of ECG events that are not wanted.

Additional parameters displayed during testing may include the following:

- Treadmill speed and grade, or watts if an ergometer is being used
- Protocol name
- Stage hold (if selected)
- Heart rate (HR)/% of target HR, and target watts if ergometer is being used
- ST level in mm or  $\mu V$  and ST slope in mV
- 3.1.12.3 • BP and SpO2 with time last acquired (optional)
- Stage and Stage Time
- Patient name
- Patient ID number
- Total exercise time
- METs and/or Double Product, and/or ST Index
- One average complex for each of the 12 leads superimposed for comparing current data to reference data
- User-defined, expanded average complex superimposed for comparing current data to reference data
- Running trends of METs with HR, BP systolic and diastolic values, and ST level

During the Recovery phase, Q-Stress will display both Patient Data and Conclusions selections that allow data entry for the final report. At the end of the Recovery phase, the Report Manager will display a summary page that allows definition and creation of the final report.

The final report may consist of the following sections that are user enabled or disabled:

- Patient Information
- Exam Summary
- Trends for Heart Rate, BP, workload, ST level, and ST slope
- Worst Case Average
- Periodic Averages
- Peak Averages
- ECG Prints
  - Automatic 12-lead ECGs per protocol
  - Peak exercise ECG
  - Arrhythmia events
  - User added 12-lead ECG Events (Supine, Standing, Symptoms, Perceived Exertion, etc.)

3.1.11.2

#### 40 Hz Filter

The 40 Hz filter is a plot filter, meaning that it affects only the plotted/printed information similar to a 40 Hz filter on an electrocardiograph. The default state of this filter (enabled or disabled) is determined by the selected profile. When the filter is on, 40 Hz appears in the lower border of the real-time ECG display. This filter setting can be toggled between 150 Hz and 40 Hz as needed during a stress test as well as during post-exam review.

It is recommended that the 40 Hz filter is disabled for most exams. This results in a standard ECG filter at 150 Hz. The 40 Hz filter may be used when poor ECG quality causes interpretation difficulties after good electrode to skin preparation has been performed.

**WARNING:** When the 40 Hz filter is used, the frequency response requirement for diagnostic ECG equipment cannot be met. The 40 Hz filter significantly reduces high-frequency components of the ECG and pacemaker spike amplitudes, and is recommended only if high-frequency noise cannot be reduced by proper procedures.

#### Beat Identification and Accuracy

At the start of a stress exam, Q-Stress initially acquires a dominant QRS complex for each lead to build the first QRS template for all 12 ECG leads. The average QRS waveform for each of the 12 leads is updated after every beat.

If the dominant QRS morphology changes, it is automatically detected and the new morphology is "learned" as the new dominant beat morphology. This event is labeled as DRC (Dominant Rhythm Change) on the displayed trends.

#### Resting ECG Interpretation



3.1.14



The Q-Stress allows for the acquisition and printing of a baseline supine 12-lead resting ECG with Mortara Instrument's VERITAS™ resting ECG interpretation algorithm. A Mason-Likar labeled resting ECG with interpretation may also be performed with torso located leads used for comparison throughout the exam.

For additional information on this algorithm, please refer to the *Physician's Guide to VERITAS with Adult and Pediatric Resting ECG Interpretation*.

D

Panel Features	Description
<p><b>Zoomed ST</b> 3.1.16</p> 	<p>One expanded average complex with superimposition of current and reference data. ST level in mm or <math>\mu</math>V and ST slope measurement in mV are also displayed. The zoomed ECG can be shown / hidden on display. The same selected ST lead is also shown in the ST change trend.</p> <p>Tick marks are displayed on the QRS complex to show isoelectric, J-point, and ST measurement points.</p> <p><i>NOTE: A click on the Zoomed ECG will open a menu allowing selection of a different lead, different reference, addition of ST baseline markers, ability to move the isoelectric and J-point for any needed correction, and a Reload of the dominant beat morphology.</i></p>
<p><b>ST Profile</b> 3.1.17</p> 	<p>Displays ST value for the on-screen average in a graph format. When Exercise begins, Q-Stress will acquire and analyze the incoming ECG data to develop the current ST level at the beginning of the procedure. The graph displays the current ST levels in black and the reference levels in green.</p> <p><i>NOTE: A click on the ST Profile will open a menu allowing change between Last Rhythm Event, ST Profile, and Profile and Event.</i></p>

Panel Features	Description
<p><b>Real-time ECG</b></p> 	<p><b>3.1.12.2</b></p> <p>3 leads, 6 leads, 12 leads (12x1), or 12 leads (6x2) are displayed in real time with their respective lead labels.</p> <p>This panel will increase/decrease in total ECG seconds displayed according to other shown panels.</p> <p>This panel is always displayed.</p> <p><i>NOTE: A click on the real-time ECG will open a menu allowing selection of displayed lead layout, Gain, Speed, and Filter change for the current patient.</i></p>
<p><b>Context View</b> 3.1.18</p>  <p>To Delete a stored ECG event:</p> <p>To Add a new ECG event:</p>	<p>Displays a miniature single ECG lead when enabled. The phase or stage together with the time is displayed in blue text.</p> <p>Choose the displayed lead by a click on the lead label when a drop-down list appears.</p> <p>A scroll bar is used to move back and forward in time from pre-exercise through the recovery phase. When scrolling back in time, the display will return to the current time after 60-second of no activity.</p> <p>Stored 12-lead ECGs are highlighted and labeled in a blue box that can be edited with a click. The blue box changes to red with an x in the upper right corner. Click on the x to delete the stored ECG.</p> <p>Double-click the ECG to add an ECG event and choose a label from the Event drop-down list or enter free text. The Bookmark label allows a quick selection that can be edited later.</p>

GOLD-STANDARD MEDICAL TREADMILLS THAT ARE AMONG THE BEST SELLING IN THE WORLD


# Quinton TM55 and TM65 Treadmills

2.2, 3.2

## PRIMARY BENEFITS:

- **Ease-of-use:** The treadmill is easy to operate and whisper-quiet at all speeds.
- **Control:** The Ultradrive system provides exceptionally smooth power and accuracy even at the slowest speeds and the heaviest weights.
- **Integration:** TM55 and TM65 treadmills are compatible with Quinton Q-Stress and Montara X-Scribe stress testing systems.

3.2.8

## THE TM55 AND TM65 PROVIDE A COMFORTABLE AND QUIET TREADMILL EXPERIENCE

Clinicians around the world rely on our treadmills. Your patients will appreciate the solid and secure feeling the TM55 and TM65 provide, while you will appreciate how easy they are to use.

- Our easy-to-read cardiac stress systems make it simple for you to start, control, and stop the treadmill.
- Clear displays let you quickly find test duration, speed, and grade.
- The stress system allows you to bring your patients to a comfortable stop, gradually reducing speed and lowering the grade to zero.
- Optional rapid deceleration capability is available for the times you need to slow and stop the treadmill quickly, such as stress echo peak exercise imaging.

## THE ULTRADRIVE SYSTEM IS THE POWER BEHIND THE TM55 AND TM65 TREADMILLS

3.2.2 Ultradrive is a digitally controlled AC drive system that provides exceptional performance, accuracy and reliability.

- Treadmills deliver high torque for accurate speed control even at lower speeds.
- Patients enjoy a natural walking/running sensation with no lag when the foot strikes the belt.

Quinton treadmills have been the industry's gold standard for more than three decades and are among the best selling in the world.

QUINTON


2,3, 3.3

Automated Blood Pressure Monitor  
for Cardiac Stress and Exercise Testing

  
**SunTech®**

Accurate, Motion Tolerant BP Monitoring for Your Stress Lab

+ Reliable, automated blood pressure,  
so you can focus on your patient

+ Dual blood pressure measurement  
modes, Exercise & Non-Exercise

SunTech Medical's Tango® M2 automated NIBP monitor uses DKA technology specifically engineered for the cardiac stress testing environment. With advanced proprietary algorithms SunTech Medical stands alone in the market of motion tolerant BP measurement. The Tango M2 allows for the focus to be on the patient by automatically taking accurate and reliable BP measurements rather than spending valuable time trying to take measurements manually. Tango M2 can add automated BP and SpO<sub>2</sub> to your stress ECG system, creating a complete cardiac diagnostic testing center.

#### + Features & Benefits + Description

	Seamless Integration:	Automated communication with your stress system reduces the risk of transcription errors.
	Non-Exercise BP Mode:	Allows BP measurement during patient set-up and recovery without an ECG signal.
	Stat Mode:	Rapidly repeated automatic BP measurements for time sensitive and emergent situations.
3.3.5	Color LCD	Improved usability with a new 7" color LCD display.
3.3.6	Verify Measurements:	See & hear the Korotkoff sounds using the onscreen display and the included headphones.
	Data Retrieval:	Easier troubleshooting with 300 BP reading history and USB capabilities for measurement retrieval.
	Field Upgrades:	USB port allows for field upgrades, making sure end-users always have the current software.
	E-Library:	Onboard E-Library, including step-by-step interface notes and tutorials.

#### + Options

Single Patient Use (SPU) kits for increased infection control  
Pulse oximetry (SpO<sub>2</sub>)  
Internal ECG

3.3.2

Orbit-K BP cuff

SpO<sub>2</sub>

#### + Specifications + Description

3.3.1	BP Measurement:	Auscultatory R-wave gating using K-sound analysis, for all static & active phases of stress testing. Oscillometric using pneumatic pressure for static measurements only.
3.3.3	Measurement Range:	Pressure (DKA Mode) Systolic: 40 - 270 mmHg, Diastolic: 20 - 160 mmHg (OSC Mode) Systolic: 40 - 260 mmHg, Diastolic: 20 - 160 mmHg; Heart Rate: 40 - 200 bpm
3.3.4	Interfaces:	Integrates with all popular stress ECG systems using RS-232, BNC, ECG cable only for internal ECG & USB connections.
	ECG Source:	Primary - From integrated stress ECG system or other external source Secondary - Internal ECG option using V2, V6, RI
	Power:	Input - 100-240 VAC @ 1.5A, 50-60 Hz. Output: +9 VDC @ 5A IEC 320 type input connector Classification - Class I, continuous
	BP Sampling Intervals:	From integrated stress ECG system or other external source, or 1 - 20 minute intervals.
	Dimensions:	24.0 cm x 17.4 cm x 11.5 cm (9.5" x 6.9" x 4.5")
	Weight:	1.68 Kg (3.725 lb) (59.6 oz)
	Warranty:	2 year standard warranty on monitor.
	Accuracy:	Equivalent to a trained observer using a cuff/stethoscope auscultation method per ANSI/AAMI/ISO 81060-2
	Standards:	IEC 60601-1:2005, IEC 60601-1-2:2007 (EMC), IEC 80601-2-30:2009, ISO 80601-2-61:2011, ISO 10993-1:2009, ISO 10993-5:2009, ISO 10993-10:2010, FDA 21CFR801.5, MDD, WEEE

**SunTech Medical**  
Brilliant Blood Pressure Solutions™

SunTech Medical, Inc.  
507 Airport Boulevard, Suite 117  
Montville, NC 27560-8200 USA  
Tel: + 1 919.654.2300  
1.800.421.8626  
Fax: + 1 919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd.  
Oakfield Estate  
Eynsham, Oxfordshire  
England  
Tel: + 44 (0) 1865.884.234  
Fax: + 44 (0) 1865.884.235

www.SunTechMed.com

SunTech Medical, Ltd.  
2/F of Building A, Jinxiongda Technology Park  
Guanlan, Bao'an District Shenzhen, 518110, PRC  
Tel: + 86-755-2958 8810  
+ 86-755-29588986 (Sales)  
+ 86-755-29588665 (Service)  
Fax: + 86-755-2958 8829

www.CN.SunTechMed.com

part #82-0126-00 R1

Hillrom™

V  
VO D.3

Hillrom Services Ptd Ltd  
1 Yishun Avenue 7  
Singapore 768923

Date: 22 November 2019

To:  
SaintMed Co.,Ltd.  
222,222/1-2 West Inn Complex Bldg.,  
10th Floor, Moo 1, Rattanaibeth Rd.,  
Bangrakyai, Bangbuathong,  
Nonthaburi 11110 Thailand

Mortara Product Spare Part Availability

Dear Valued Partner,

We hereby confirm that the following product spare part are available upon 5 years from End-of-life date.

- Q-stress
- X-Scribe
- Treadmill TM55
- Treadmill TMX428
- H-Scribe
- H3+

Thank you.

Name: Yong Yuen Shen  
Head of Service, APAC  
Date: 22 November 2019

Hillrom Services Pte. Ltd. (U.S. / GST Reg. No.: 2637235824)  
1 Yishun Avenue 7  
Singapore 768923  
hillrom.com  
Tel: +65 6199 7350  
Fax: +65 6199 7351

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